



Freedom Designs, Inc.

# Freedom SP3—Mini Series

2241 N. Madera Rd. Simi Valley, CA 93065 Phone: (800) 331-8551 Fax (888) 582-1509 www.freedomdesigns.com

|                 |                        |                 |                                     |
|-----------------|------------------------|-----------------|-------------------------------------|
| <b>Bill to:</b> | Date: _____ PO#: _____ | <b>Ship to:</b> | Company: _____                      |
|                 | Company: _____         |                 | Address: _____                      |
|                 | Buyer: _____           |                 | City: _____ State: _____ Zip: _____ |
|                 | Marked For: _____      |                 | Phone: _____                        |


| Description | Part No. | Retail | Description | Part No. | Retail |
|-------------|----------|--------|-------------|----------|--------|
|-------------|----------|--------|-------------|----------|--------|

**Frame Type**


**Freedom SP3 Mini Folding** SPM-30270 **\$2315**  
 (E1238) or (K0005)  
 HCPCS codes for frames 10"W—18"W

With double seat rail for upholstery & seating.  
 Without double seat rail for seating only. \*

Folding Frame  
With Seat Sling




Folding Frame  
For Seating Only




**Freedom SP3 Mini Rigid** SPM-30270R **\$2315**  
 (E1237) or (K0005)  
 HCPCS codes for frames 10"W—18"W

With double seat rail for upholstery & seating.  
 Without double seat rail for seating only. \*

Rigid Frame  
With Seat Sling



Rigid Frame  
For Seating Only



\*Seat to floor height - 1 5/8" less  
 \*Lighter weight

**Transport Option**

**Transport Package** SPM-99998 **\$150**  
 Available with FDI WC19 approved seating  
 or on frame only.


**No Transport Package** SPM-99999-1 **n/c**

**Configuration Set-up**


**Standard—Rear Wheel Drive** SPM-STD **n/c**

**Reverse Configuration—Front Wheel Drive** SPM-REV **n/c**

Standard  
Configuration




Reverse  
Configuration



**Frame Width**

Specify Frame Width:

|  |           |       |
|--|-----------|-------|
| <input type="radio"/> <b>10" Wide</b> —Available without sling only. | SPM-30250 | n/c   |
| <input type="radio"/> <b>11" Wide</b>                                | SPM-30210 | n/c   |
| <input type="radio"/> <b>12" Wide</b>                                | SPM-30220 | n/c   |
| <input type="radio"/> <b>13" Wide</b>                                | SPM-30230 | n/c   |
| <input type="radio"/> <b>14" Wide</b>                                | SPM-30240 | n/c   |
| <input type="radio"/> <b>15" Wide</b>                                | SPM-30290 | n/c   |
| <input type="radio"/> <b>16" Wide</b>                                | SPM-30260 | n/c   |
| <input type="radio"/> <b>17" Wide</b>                                | SPM-30330 | \$ 75 |
| <input type="radio"/> <b>18" Wide</b>                                | SPM-30280 | \$ 75 |



← Width →

**Frame Length**

Specify Frame Length:

|   |           |      |
|---|-----------|------|
| <input type="radio"/> <b>Small (9" - 12")</b> | SPM-30509 | n/c  |
| <input type="radio"/> <b>9" Deep</b>          | SPM- 9    | n/c  |
| <input type="radio"/> <b>10" Deep</b>         | SPM-10    | n/c  |
| <input type="radio"/> <b>11" Deep</b>         | SPM-11    | n/c  |
| <input type="radio"/> <b>12" Deep</b>         | SPM-12    | n/c  |
| <input type="radio"/> <b>Long (11" - 14")</b> | SPM-30510 | \$75 |
| <input type="radio"/> <b>11" Deep</b>         | SPM-11    | n/c  |
| <input type="radio"/> <b>12" Deep</b>         | SPM-12    | n/c  |
| <input type="radio"/> <b>13" Deep</b>         | SPM-13    | n/c  |
| <input type="radio"/> <b>14" Deep</b>         | SPM-14    | n/c  |

**Sling Seat**

|  |           |               |
|--|-----------|---------------|
| <input type="radio"/> <b>Omit Seat Sling</b> | SPM-39004 | <b>(\$15)</b> |
| <input type="radio"/> <b>Seat Sling</b>      | SPM-39011 | n/c           |

Specify Sling Length:  
 Front of seat rail to back upholstery.

**Seat Sling Length:**

|                                       |           |     |
|---------------------------------------|-----------|-----|
| <input type="radio"/> <b>9" Deep</b>  | SPM-30409 | n/c |
| <input type="radio"/> <b>10" Deep</b> | SPM-30500 | n/c |
| <input type="radio"/> <b>11" Deep</b> | SPM-30501 | n/c |
| <input type="radio"/> <b>12" Deep</b> | SPM-30502 | n/c |
| <input type="radio"/> <b>13" Deep</b> | SPM-30503 | n/c |
| <input type="radio"/> <b>14" Deep</b> | SPM-30504 | n/c |

**Frame Colors—Powder Coated**

|                                     |                                      |
|-------------------------------------|--------------------------------------|
| <input type="radio"/> Black         | <input type="radio"/> Pacific Blue   |
| <input type="radio"/> Blizzard Blue | <input type="radio"/> Passion Purple |
| <input type="radio"/> Blue Blast    | <input type="radio"/> Pink Parfait   |
| <input type="radio"/> Blue Stardust | <input type="radio"/> Really Red     |
| <input type="radio"/> Cosmic Purple | <input type="radio"/> Red Sunset     |
| <input type="radio"/> Extreme Green | <input type="radio"/> Sparkle Green  |
| <input type="radio"/> Mystic Teal   | <input type="radio"/> Twilight       |
| <input type="radio"/> Octane Orange | <input type="radio"/> Yellow         |

**Cushions**

|   |           |              |
|---|-----------|--------------|
| <input type="radio"/> <b>1" Cushion</b> | SPM-39003 | <b>\$ 85</b> |
| <input type="radio"/> <b>2" Cushion</b> | SPM-39007 | <b>\$ 95</b> |

Includes 4# foam & black knit cover.

Date: \_\_\_\_\_ PO#: \_\_\_\_\_

Marked For: \_\_\_\_\_

| Description | Part No. | Retail | Description | Part No. | Retail |
|-------------|----------|--------|-------------|----------|--------|
|-------------|----------|--------|-------------|----------|--------|

**Back Canes / Push Handles  
Hand Grip Bend Options**

- Traditional Bend Push Handles** SPM-30518 n/c
- Ergo Bend Push Handles** SPM-30519 n/c



Traditional Push Handle

Ergo Push Handle

Specify height from rear seat rail to top of bend on cane:

- Short (19" High)** SPM-30515 n/c
- Medium (23" High)** SPM-30514 n/c
- Tall (26" High)** SPM-30516 n/c

- Push Handle Extension - Quick Release** SPM-30517 \$100
- One Piece Push Handle Extension - Quick Release** SPM-30526 \$125
- One Piece Push Button Angle Adjustable Stroller Handle** SPM-30527 \$175

Available on chairs 14"W and up.  
Note: Chair will not fold with this stroller handle.

Back canes can achieve an additional height of 6.5" when quick release push handle extensions are ordered.

**Sling Backs**

- Omit Sling Back** SPM-39005 (\$15)
- Sling Back\*** SPM-39015 n/c
- Tension Adjustable Back\*** SPM-39008 \$231

Specify Sling Back Height

(Top of seat to top of back upholstery)

- 10" High** SPM-30910 n/c
- 11" High** SPM-30911 n/c
- 12" High** SPM-30912 n/c
- 13" High** SPM-30913 n/c
- 14" High** SPM-30914 n/c
- 15" High** SPM-30915 n/c
- 16" High** SPM-30916 n/c
- 17" High** SPM-30917 n/c
- 18" High** SPM-30918 n/c


\*Back upholstery height must be chosen when ordering sling or tension adjustable back upholstery

**Swing Away Hangers**

- Omit Hangers & Footplates** SPM-33000 (\$30)
- Spring Loaded Pop-Up** (pictured below) SPM-30272 n/c
- Latch Style Swing Away** SPM-30278 n/c  
*Not available in reverse configuration.*

Specify:

- 90 Degree Hangers\*** SPM-33334 n/c
  - 75 Degree Hangers\*** SPM-33335 n/c
- \*See Configuration Chart



90° Latch Style Hanger

75° Pop-Up Style Hanger

**Footrest Options**



Angle Adjustable Flip Up Footplate

Removable Footboard

- Omit Footplates Only** SPM-36404 (\$15)
- Composite Footplate - Flip Up** SPM-17070 n/c
- Kids Footplate (K0040) Flip Up Footplate—Adjustable Angle** SPM-36020 \$150  
For chairs 14" wide and under.
- Adult Footplate 15"W & Over (K0040) Flip Up Footplate—Adjustable Angle w/ Heel Loop** SPM-36021 \$150  
For chairs 15" wide and over.
- Omit Heel Loop** SPM-36023 n/c
- Heel Loop (E0951)** SPM-36022 \$ 50
- Footboard—Removable, Angle Adj.** SPM-36029 \$150
- Hanger Extensions - .5—4.5" Long** SPM-36402 \$ 55
- Hanger Extensions - 1.5—6.5" Long** SPM-36403 \$ 65  
1.5" - 6.5" long only available with 75 degree hangers.

| Configuration                     | Knee to Heel Range (from seat rail) |
|-----------------------------------|-------------------------------------|
| Hanger with High Mount Clamps     | 4" - 9"                             |
| Hanger Only                       | 6" - 9"                             |
| Hanger with .5" - 4.5" Extension  | to 13.5"                            |
| Hanger with 1.5" - 6.5" Extension | to 15.5"                            |

- High Mount Clamp (K0037)** SPM-36407 \$100  
Not available on 75 degree hangers.
- Knee to Heel Measurement** \_\_\_\_\_ Inches

Date: \_\_\_\_\_ PO#: \_\_\_\_\_

Marked For: \_\_\_\_\_

| Description | Part No. | Retail | Description | Part No. | Retail |
|-------------|----------|--------|-------------|----------|--------|
|-------------|----------|--------|-------------|----------|--------|

**Casters**

- Omit Caster Wheels** SPM-34333 n/c
- 3" Poly Casters** SPM-34888 \$ 55
- 3" Light Up Casters** SPM-34111 \$ 55  
3" casters in combination with 20" wheels are not available.  
3" casters recommended for low frame configurations.
- 4" Poly Casters** SPM-40000 \$ 65
- 4" Light Up Casters** SPM-41000 \$ 99
- 5" Poly Casters** SPM-34555 n/c
- 5" x 1.5" Poly Casters** SPM-34222 \$ 85
- 6" Poly Casters** SPM-34666 n/c
- 6" x 1.5" Poly Casters** SPM-34999 \$ 85
- 6" x 2" Poly Casters with Tread** SPM-34444 \$ 90
- 6" Pneumatic Casters** SPM-34667 n/c
- Frog Legs (E1015)** SPM-35000 \$325  
Not available with 6" x 2" casters.
- Caster Stem Bolt—.75" Longer** SPM-34778 n/c
- Caster Stem Bolt—1.5" Longer** SPM-34779 n/c



**Anti-Tip Tubes**

- Anti-Tip Tubes—Fixed (E0971)** SPM-31600 \$ 95
- Anti-Tip Tubes—Swing Away (E0971)** SPM-31602 \$ 95  
Can not be used with reverse configuration.

**Wheels & Wheel Accessories**

Patented double lock axle receiver included at no charge.

- Omit Wheels & Axles** SPM-31524 (\$55)

Mag Wheel



**Mag Wheels**

Spoke Wheel



- 12" Pneumatic Mag** SPM-31525 n/c
- 12" Polyurethane Mag** SPM-31526 \$ 35
- 12" Semi Pneumatic Mag** SPM-31523 \$165  
Handrims not available with 12" wheels.
- 16" Pneumatic Mag** SPM-31530 \$ 40
- 18" Pneumatic Mag** SPM-31533 \$ 40
- 20" Pneumatic Mag** SPM-31527 n/c  
Available on long frame only.  
Aluminum handrims are included on 20" mag wheels.  
To omit handrims, please check below

**Spoke Wheels**

- 20" Pneumatic Spoke** SPM-31540 \$200  
Available on long frame only.  
Aluminum handrims are included on 20" spoke wheels.  
Plastic coated handrims are not available with spoke wheels.  
Spoke wheels are not available with hub locks.

**Wheel Accessories**

- Omit Handrims** SPM-31006 (\$20)  
On 20" mag or spoke wheels.
- Aluminum Handrims 16" & 18" (pr.)** SPM-31007 \$166
- Plastic Coated Handrims (pr.)** SPM-31005 \$117  
Plastic coated handrims are not available with spoke wheels.
- Natural Fit Handrims**  
Please contact your customer service representative for details.
- Airless Inserts (pr.) (E2213)** SPM-31000 \$ 75

**Wheel Locks**

Choose type of lock - either standard or hub locks.

Specify lever type:

- Push to Lock\*** SPM-31501 n/c
- Pull to Lock\*** SPM-31506 n/c
- Rear Mount Push to Lock** SPM-31507 n/c
- Rear Mount Pull to Lock** SPM-31508 n/c
- Wheel Lock Handle Ext. (E0961)** SPM-31801 \$ 60



\* If both wheel locks and hub locks are ordered, an "off chair" charge will be applied.

**OR**

**Hub Locks**

(Not available with spoke wheels.)

Specify lever location on frame:

Determines cable length and location of engagement lever for client access.

- Front Frame** SPM-31509 n/c
- Middle Frame** SPM-31510 n/c
- Rear Frame** SPM-31511 n/c
- Back Canes** SPM-31512 n/c

**Armrests**

- Omit Flip Back Armrest** SPM-38209 (\$20)

**Back Cane Mount**

- Flip Back Armrests—13" Long** SPM-38204 n/c
- Flip Back Armrests—16" Long** SPM-38202 n/c
- Locking Flip Back Armrests—13" Long (E0973)** SPM-38201 \$235
- Locking Flip Back Armrests—16" Long (E0973)** SPM-38200 \$235
- Flip Back Armrests—Telescoping** GZ-38208 \$100
- Locking Flip Back Armrests—Telescoping** GZ-38207 \$260

**Frame Mount**

- Height Adjustable Padded Armrest (E0973)** SPM-38203 \$235



Flip Back Armrest



Height Adjustable Padded Armrest



Locking Flip Back Armrest


- Desk Length Flat Armrest Pad** SPM-38213 n/c
- Full Length Flat Armrest Pad** SPM-38214 n/c  
Only available on 16" armrests.
- Seat Rail to Top of Armrest** \_\_\_\_\_ Inches

Date: \_\_\_\_\_ PO#: \_\_\_\_\_

Marked For: \_\_\_\_\_

| Description | Part No. | Retail | Description | Part No. | Retail |
|-------------|----------|--------|-------------|----------|--------|
|-------------|----------|--------|-------------|----------|--------|

**Seating Options**



- Rigid Seat Base SPM-39006 \$100
- Removable Seat Base SPM-39012 \$150
- Seat Rail to Floor Height \_\_\_\_\_ Inches


**Accessories**

- Medical Necessities Bag (black only) SPM-23207 \$67
- Back Pack (black only) SPM-TFBP \$55

**Canopy**

- Canopy Fabric with Frame\* SPM-38750 \$150
- Canopy with Visor and Frame\* SPM-38756 \$175

\*Canopies w/transport option must be removed during transport.



**Notes:**

**\*Positioning Belts (E0978)**

- Streamline - Two Point Padded**
- Center Pull - Metal Push Button Buckle SB-3332 \$ 88
  - Center Pull - Metal Push Button with 1" Buckle & .75" Webbing SB-3332-1 \$ 88
  - Center Pull - EZ Release Push Button Buckle SB-3334 \$ 88
- Streamline - Non-Padded**
- Center Pull - Metal Push Button Buckle SB-3410 \$ 29
- Hipster - Two Point Padded**
- Center Pull - Metal Push Button Buckle HIP-3002 \$109
  - Center Pull - EZ Release Push Button Buckle HIP-3004 \$109
- Hipster - Non-Padded**
- Center Pull - Metal Push Button HIP-3021 \$ 50
  - Center Pull - Metal Push Button with 1" Buckle & .75" Webbing HIP-3021-1 \$ 50

\*Many other belts are available—those shown above include some of our most popular. Ask your customer service representative about additional styles.

**Wheel & Caster Matrix**

| Floor to Front of Seat Rail | Caster Size | Wheel Size | Arc Over Rail | Fork Size |
|-----------------------------|-------------|------------|---------------|-----------|
| 10.5"                       | 3"          | 16"        | 4.5"          | Short     |
| 11.5"                       | 3"          | 18"        | 6.5"          | Short     |
| 12.5"                       | 5"          | 18"        | 5.5"          | Long      |
| 12.5"                       | 5"          | 20"        | 7.5"          | Long      |
| 13"                         | 5"          | 18"        | 5"            | Long      |
| 13"                         | 5"          | 20"        | 7"            | Long      |
| 13.5"                       | 5"          | 18"        | 4.5"          | Long      |
| 13.5"                       | 5"          | 20"        | 6.5"          | Long      |

**Maximum Rider Weight—200 lbs.  
Frame weight is approximately 23 lbs.**