



# Freedom CGX™

## Center of Gravity eXperience



- A full range of adjustment from  $-10^{\circ}$  anterior to  $+50^{\circ}$  posterior
- Seat widths and depths of 10" - 22"
- Foot-activated tilt

**Now Available!**  
**250 lb weight capacity on**  
**10" - 18" widths and depths**

More than mobility products,  
we make a difference  
one client at a time.



### **Aesthetics**

Sleek refined appearance.

### **Durability**

A strong rigid frame built around a simple and robust tilt design.

### **Dependability**

Utilizes no tilt cables or mechanical locks.

### **Versatility**


Adjusts 4" in width and 4" in depth accommodating client needs.

### **Adaptability**

Built-in growth.



# Freedom CGX™



**Back posts fold easily for convenient transportation**



**Foot activated tilt**



**New latch style hanger**

**-10° anterior tilt**  
\*Not available on Extra Wide



## Specifications

<b>Frame Style</b> Rigid, Tilt-in-Space	<b>Wheel Sizes</b> 12", 16", 18", 20", 22", 24"
<b>Seat Widths</b> Narrow 10" - 14" Medium 12" - 16" Wide 14" - 18" Extra Wide 18" - 22"	<b>Hanger Options</b> Pop up removable: 75° and 90° Latch Style Swing Away: 75° and 90°
<b>Seat Depth</b> Short 10" - 14" Medium 12" - 16" Long 14" - 18" Extra Wide 18" - 22"	<b>Adjustable AngleBack</b> 80°, 90°, 100°, 110°, 120°
<b>Back Height</b> Short: 21" Medium: 24" Tall: 27"	<b>Tilt Range</b> -10° to +50° 0° to +48° X-Wide
<b>Rider Weight Limit</b> 250 lbs - N, M, W 300 lbs - Extra Wide	<b>Frame Weight</b> 38 lbs
<b>Caster Sizes</b> 3", 4", 5", 6", 7"	<b>• Extra Wide includes lift assist for stable and nearly effortless tilt for heavier clients</b>

2241 N. Madera Rd. Simi Valley, CA 93065  
Phone: 800-331-8551 • Fax: 888-582-1509  
Email: CustomerService@freedomdesigns.com