



Bill to:	Date: _____ PO#: _____	Ship to:	Company: _____
	Company: _____		Address: _____
	Buyer: _____		City: _____ State: _____ Zip: _____
	Marked For: _____		Phone: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Frame Type		
<input type="radio"/> TriPod LT Pediatric (Rigid) (E1237)	LPTP-30270	\$2315
Front Frame Configuration		
<input type="radio"/> 90 Degree Straight	LPTP-30330	n/c



Frame Length			
<u>Specify Frame Length:</u>			
<input type="radio"/> Small (11" - 13")	LPTP-30507	n/c	
<input type="radio"/> 11" Deep	LPTP-11	n/c	
<input type="radio"/> 12" Deep	LPTP-12	n/c	
<input type="radio"/> 13" Deep	LPTP-13	n/c	
Shorter seat depth may be achieved through seating system. Seat depth indicates back cane placement.			

Transport Option		
<input type="radio"/> Transport Package	LPTP-99999	\$150
Available with FDI WC19 approved seating or on frame only.		
<input type="radio"/> No Transport Package	LPTP-99999-1	n/c

Seating Options			
<input type="radio"/> Omit Seat Pan or Sling	LPTP-39011	n/c	
<input type="radio"/> Seat Pan - Fixed	LPTP-39010	\$100	
<input type="radio"/> Seat Pan w/EZ Lite	LPTP-39013	\$150	
<input type="radio"/> Seat Sling	LPTP-39012	n/c	

Frame & Axle Colors		
<input type="radio"/> Black	<input type="radio"/> Pacific Blue	
<input type="radio"/> Blizzard Blue	<input type="radio"/> Passion Purple	
<input type="radio"/> Blue Blast	<input type="radio"/> Pink Parfait	
<input type="radio"/> Blue Stardust	<input type="radio"/> Really Red	
<input type="radio"/> Cosmic Purple	<input type="radio"/> Red Sunset	
<input type="radio"/> Extreme Green	<input type="radio"/> Sparkle Green	
<input type="radio"/> Mystic Teal	<input type="radio"/> Twilight	
<input type="radio"/> Octane Orange	<input type="radio"/> Yellow	
Axle Colors		
<input type="radio"/> Black	<input type="radio"/> Yellow	<input type="radio"/> Same as Frame

<u>Specify Sling Length:</u> Front of seat rail to back upholstery.			
<input type="radio"/> Seat Sling Length:			
<input type="radio"/> 11" Deep	LPTP-30501	n/c	
<input type="radio"/> 12" Deep	LPTP-30502	n/c	
<input type="radio"/> 13" Deep	LPTP-30503	n/c	
<input type="radio"/> 14" Deep	LPTP-30504	n/c	

Frame Widths		
<u>Specify Frame Width:</u>		
<input type="radio"/> 10" Wide	LPTP-30200	n/c
<input type="radio"/> 12" Wide	LPTP-30220	n/c
<input type="radio"/> 14" Wide	LPTP-30240	n/c
<input type="radio"/> 16" Wide	LPTP-30260	n/c
<input type="radio"/> 18" Wide	LPTP-30280	n/c

Cushions			
<input type="radio"/> 1" Cushion	LPTP-39003	\$ 85	
<input type="radio"/> 2" Cushion	LPTP-39007	\$ 95	
Includes 4# foam & black knit cover.			

Date: _____ PO#: _____

Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

**Back Canes / Push Handles
Hand Grip Bend Options**

- Straight** - No push handles
With fold over sling back LPTP-30600 n/c
- Specify height from rear seat rail to top of straight cane:
- 10" High** LPTP-30610 n/c
- 11" High** LPTP-30611 n/c
- 12" High** LPTP-30612 n/c
- 13" High** LPTP-30613 n/c
- 14" High** LPTP-30614 n/c
- 15" High** LPTP-30615 n/c
- 3/4" Telescoping Push Handles** LPTP-30616 **\$175**
(10" - 15") These are push handles only.
Back upholstery does not attach to telescoping handles.
For use with straight back posts only.
- Traditional Bend Push Handles** LPTP-30518 n/c
- Ergo Bend Push Handles*** LPTP-30519 n/c



Traditional Push Handle



Ergo Push Handle

- Specify height from rear seat rail to top of bend on cane:
 - Short (19" High)** LPTP-30515 n/c
 - Medium (23" High)** LPTP-30514 n/c
 - Tall (26" High)** LPTP-30516 n/c
 - Push Handle Extension-Quick Release** LPTP-30517 **\$100**
 - One Piece Push Handle Extension - Quick Release** LPTP-30526 **\$125**
 - One Piece Push Button Angle Adjustable Stroller Handle** LPTP-30527 **\$175**
Available on chairs 14"W and up.
- Back canes can achieve an additional height of 6.5" when either quick release push handle extensions are ordered.

Back Upholstery Height

- Omit Sling Back** LPTP-39005 n/c
- Sling Back*** LPTP-39015 n/c
- Tension Adjustable Back*** LPTP-30008 **\$231**
- Specify Sling Back Height:
(Measured from top of seat to top of back upholstery)
- 10" High** LPTP-30910 n/c
- 11" High** LPTP-30911 n/c
- 12" High** LPTP-30912 n/c
- 13" High** LPTP-30913 n/c
- 14" High** LPTP-30914 n/c
- 15" High** LPTP-30915 n/c
- 16" High** LPTP-30916 n/c
- 17" High** LPTP-30917 n/c
- 18" High** LPTP-30918 n/c

*Back upholstery height must be chosen when ordering sling or tension adjustable back upholstery.

Footrest Options

- Angle Adjustable, Flip-Up Footboard** LPTP-36028 **\$125**
(K0040)
Adjustable leg length from 3.5" - 12"
Flip up footboard available with straight (non taper) frames only.
- Individual Flip Up/Angle Adjustable Footplates** (K0040) LPTP-36030 **\$105**
- Composite Individual Flip Up Footplates** LPTP-36031 n/c
Adjustable leg length from 4.5" - 6.5".
All measurements are without a cushion.



Angle Adjustable Flip Up Footplate



Flip Up Footboard

Casters



- Omit Casters Wheels** LPTP-34333 n/c
- 3" Poly Casters** LPTP-34888 **\$ 55**
- 3" Light Up Poly Casters** LPTP-34111 **\$ 55**
- 4" Poly Casters** LPTP-40000 n/c
- 4" Light Up Poly Casters** LPTP-41000 **\$ 99**

Center Wheel Camber Angle

- Specify:
- 0 Degrees** LPTP-31010 n/c
- 3 Degrees** LPTP-31011 n/c
- 6 Degrees** LPTP-31012 n/c



Date: _____ PO#: _____

Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Wheels & Wheel Accessories

Patented double lock axle receiver included at no charge.

- Omit Wheels & Axles** LPTP-31524 (\$55)

Mag Wheels



- 16" Pneumatic Mag** LPTP-31530 n/c
- 18" Pneumatic Mag** LPTP-31533 n/c

Wheel Accessories

- Aluminum Handrims 16" & 18" (pr.)** LPTP-31007 n/c
- Natural Fit Handrims**
Please contact your customer service representative for details.
- Airless Inserts (pr.) (E2213)** LPTP-31000 \$ 75

Wheel Locks

Choose one type of wheel lock - either standard or hub locks.

Specify lever type:



- Push to Lock*** LPTP-31501 n/c
- Pull to Lock*** LPTP-31506 n/c
- Wheel Lock Handle Ext. (E0961)** LPTP-31801 \$ 60

*If both wheel locks and hub locks are ordered, an "off-chair" charge will be applied.

OR

Hub Locks

(Not Available with spoke wheels.)

Specify Lever Location:

Determines cable length and location of engagement lever for client access.

- Front Frame** LPTP-31509 n/c
- Middle Frame** LPTP-31510 n/c
- Rear Frame** LPTP-31511 n/c
- Back Canes** LPTP-31512 n/c

Armrests

- Omit Flip Back Armrests** LPTP-38209 (\$20)
- Back Cane Mount**
- Flip Back Armrests—13" Long** LPTP-38204 n/c
 - Flip Back Armrests—16" Long** LPTP-38202 n/c
 - Locking Flip Back Armrests—13" Long (E0973)** LPTP-38201 \$235
 - Locking Flip Back Armrests—16" Long (E0973)** LPTP-38200 \$235
 - Flip Back Armrests—Telescoping** LPTP-38208 \$100
 - Locking Flip Back Armrests—Telescoping** LPTP-38207 \$260

Frame Mount

- Height Adjustable Padded Armrests (E0973)** LPTP-38203 \$235
Not available with sling or seat pan.



Flip Back Armrest



Locking Flip Back Armrest



Height Adjustable Padded Armrest

- Desk Length Flat Armrest Pad** LPTP-38213 n/c
- Full Length Flat Armrest Pad** LPTP38214 n/c
Only available on 16" armrests.
- Seat Rail to Top of Armrest** _____ Inches

Accessories

- Medical Necessities Bag** LPTP-23207 \$ 67
- Back Pack (black only)** LPTP-TFBP \$ 55
- Clothing Guards (fabric)** LPTP-CGF \$120
- Clothing Guards (Quick Release - rigid)** LPTP-CGR \$120

Date: _____ PO#: _____

Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Positioning Belts		
Streamline - Two Point Padded		
<input type="radio"/> Center Pull - Metal Push Button Buckle	SB-3332	\$ 88
<input type="radio"/> Center Pull - Metal Push Button with 1" Buckle & .75" Webbing	SB-3332-1	\$ 88
<input type="radio"/> Center Pull - EZ Release Push Button Buckle	SB-3334	\$ 88
Streamline - Non-Padded		
<input type="radio"/> Center Pull - Metal Push Button Buckle	SB-3410	\$ 29
Hipster - Two Point Padded		
<input type="radio"/> Center Pull - Metal Push Button Buckle	HIP-3002	\$109
<input type="radio"/> Center Pull - EZ Release Push Button Buckle	HIP-3004	\$109
Hipster - Non-Padded		
<input type="radio"/> Center Pull - Metal Push Button	HIP-3021	\$ 50
<input type="radio"/> Center Pull - Metal Push Button with 1" Buckle & .75" Webbing	HIP-3021-1	\$ 50
*Many other belts are available—those shown above include some of our most popular. Ask your customer service representative about additional styles.		

Wheel & Caster Matrix			
Wheel Size	Caster Size	Floor to Front of Seat Rail	Arc Over Rail
16"	3"	11"	4.5"
18"	3"	11.5"	5.5"
16"	4"	11.5"	4.5"
18"	4"	12"	5.5"
Maximum Rider Weight—200 lbs. Frame weight approximately 16.5 lbs.			

Notes:

