



Freedom[®] NXT[™]

Generation NEXT **Mini** Tilt-in-Space Mobility

Freedom: for the littlest Kids
(and their best friends too)



Folding model

NEW height adjustable
stroller push handles

Super
lightweight
yet strong
& dependable

Easy to adjust
for growing kids

Rigid model

Frame widths as narrow as 10" and frame depths as short as 8"

Good things do come in small packages!

NXT[™]Mini
Generation NEXT

What makes the Mini so mighty...

For little kids?

The NXT Mini is designed to fit the "little" kids more precisely. The seat depth adjusts 4" in ½" increments from only 8" to 12." With Freedom's "Very First" Seating System, the Mini can accommodate as little as a 5" seat depth. Foot rest heights for accurate foot positioning can start at 3" from the top of the seat frame and adjust in ½" increments to 10" with a 4" extension. The new flip back armrest is only 8.5" long to support little arms and we now have smaller foot plates for supporting little feet.



The Mini is designed with little people's safety in mind

For Moms and Dads?

The NXT Mini is packed with features to make using a wheelchair safe and simple. Freedom Designs has the only attendant hand controlled combination tilt and braking system. This system controls the speed on un-level surfaces, helps turn the wheelchair, position the seat frame in and out of tilt and locks the wheels all from the attendant's location. Height adjustable stroller push handles fit the tallest moms and the shortest dads and when the frame is folded for transport, they fold down flat so the Mini fits in tiny spaces. There is also a curb assist anti-tipper so that when navigating over small obstacles, elevating the front wheels is easy. For families on the go the Mini folds into a super lightweight compact package that is easy to lift and stow.

NXT Mini Specifications

Seat Frame Width:

Adjusts in width without parts from 10" to 12" in 1" increments

Seat Frame Depth:

Adjusts in depth without parts from 8" to 12" in ½" increments

Seat Frame Heights:

15" to 19" with various combinations of casters and wheels

Foot Plate Heights:

3" to 10" in ½" increments with various mounts or extensions

Back Post Heights:

12" to 20" & 14" to 22" with height adjustable stroller push handles

Tilt Range:

0° to 45° with ability to set max tilt angle at anything in between

Tilt Operations:

Single hand lever, dual hand lever, foot activated or combination hand activated tilt, brake and wheel lock system

Brake Operations:

Standard manual, foot activated, hand activated or combination hand activated tilt, brake and wheel lock system

Rear Wheels:

12", 16" & 18" Mag and 16" & 18" Spoke

Front Casters:

3", 4", 5", 6" & 7" in various configurations

Rider Weight Capacity:

150 lbs.

Frame Transport Weight:

19 lbs. in its folded transport configuration

Frame Colors:

21 fun colors with seating upholstery to match



Our kids are our passion.™

Freedom Designs Incorporated
2241 N. Madera Road, Simi Valley, CA 93065

phone: 800-331-8551

fax: 888-582-1509

e: CustomerService@freedomdesigns.com

www.freedomdesigns.com

Freedom Designs Incorporated (FDI) maintains a policy of continual product improvement and reserves the right to change features, specifications, and prices without notification. Check with FDI for the latest information.

Copyright © 2014 Freedom Designs, Inc.
All rights reserved.
Printed in the USA.
WBNXM 12014



Manufactured in the U.S.A.