



Freedom 2™ Freedom 2 Kids™



- Lightweight Manual Folding Wheelchair
- Can accommodate Sling Seat and Back or a Full Seating System
- Depth Adjustable Sling Back



More than mobility products,
we make a difference
one client at a time.

Aesthetics
Sleek refined
appearance.

Customizable
Sling upholstery
or custom solid
seating.


Safety
WC-19 crash
tested.

Versatility
Rigid or Folding
frame.

Stability
Engineered to
limit tipping.



Freedom 2™ Freedom 2 Kids™



Transport Option
Brackets

Freedom 2 Specifications

Frame Style Rigid, Folding	Cushions 1" or 2"
Seat Widths 11" - 18"	Back Height Short: 19" Medium: 23" Tall: 26"
Seat Depth Regular 15" - 19" Long 17" - 21"	Wheel Sizes 12", 20", 22", 24"
Sling Seat 15" - 21"	Caster Sizes 3", 4", 5", 6", 7"
Sling Back 10" - 18"	Hanger Options Pop up removable: 75° and 90° Latch Style Swing Away Elevating Legrest
Frame Weight 16.5 lbs	
Rider Weight 250 lbs	

Freedom 2 Kids Specifications

Frame Style Rigid, Folding	Cushions 1" or 2"
Seat Widths 11" - 18"	Back Height Short: 19" Medium: 23" Tall: 26"
Seat Depth Regular 10" - 14" Long 13" - 17"	Wheel Sizes 12", 20", 22", 24"
Sling Seat 10" - 17"	Caster Sizes 3", 4", 5", 6", 7"
Sling Back 10" - 18"	Hanger Options Pop up removable: 75° and 90° Latch Style Swing Away Elevating Legrest
Frame Weight 16.5 lbs	
Rider Weight 250 lbs	

Hub Lock Feature



2241 N. Madera Rd. Simi Valley, CA 93065
Phone: 800-331-8551 • Fax: 888-582-1509
Email: CustomerService@freedomdesigns.com