



Freedom CGX™

Center of Gravity eXperience



- A full range of adjustment from -10° anterior to $+50^{\circ}$ posterior
- Seat widths and depths of 10" - 22"
- Foot-activated tilt

More than mobility products,
we make a difference
one client at a time.



Aesthetics

Sleek refined appearance.

Durability

A strong rigid frame built around a simple and robust tilt design.

Dependability

Utilizes no tilt cables or mechanical locks.

Versatility


Adjusts 4" in width and 4" in depth accommodating client needs.

Adaptability

Built-in growth.



Freedom CGX™



Back posts fold easily for convenient transportation



Foot activated tilt



New latch style hanger



-10° anterior tilt
*Not available on Extra Wide

Specifications

Frame Style Rigid, Tilt-in-Space	Wheel Sizes 12", 16", 18", 20", 22", 24"
Seat Widths Narrow 10" - 14" Medium 12" - 16" Wide 14" - 18" Extra Wide 18" - 22"	Hanger Options Pop up removable: 75° and 90° Latch Style Swing Away: 75° and 90°
Seat Depth Short 10" - 14" Medium 12" - 16" Long 14" - 18" Extra Wide 18" - 22"	Adjustable AngleBack 80°, 90°, 100°, 110°, 120°
Back Height Short: 21" Medium: 24" Tall: 27"	Tilt Range -10° to +50° 0° to +48° X-Wide
Rider Weight Limit 250 lbs - N, M, W 300 lbs - Extra Wide	• Extra Wide includes lift assist for stable and nearly effortless tilt for heavier clients
Frame Weight 38 lbs	
Caster Sizes 3", 4", 5", 6", 7"	

2241 N. Madera Rd. Simi Valley, CA 93065
Phone: 800-331-8551 • Fax: 888-582-1509
Email: CustomerService@freedomdesigns.com