



Freedom Designs, Inc.

# FREEDOM 2 KIDS

2241 N. Madera Rd. Simi Valley, CA 93065 Phone: (800) 331-8551 Fax (888) 582-1509 www.freedomdesigns.com

**Bill to:**  
 Company: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_ Zip \_\_\_\_\_  
 Phone: \_\_\_\_\_

**Ship to:**  
 Company: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_ Zip \_\_\_\_\_  
 Phone: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Frame Type		
<input type="radio"/> <b>Freedom 2 Kids Depth Adjustable Rear Frame</b> (E1236) or (K0005)	FD2K-30270	<b>\$2085</b>
<u>Specify Frame Type:</u>		
<input type="radio"/> <b>Folding</b>	FD2K-30275	n/c
<input type="radio"/> <b>Rigid</b>	FD2K-30276	n/c
Transport Option		
<input type="radio"/> <b>Transport Package</b> Available with FDI WC19 approved seating or on frame only.	FD2K-99998	<b>\$150</b>
<input type="radio"/> <b>No Transport Package</b>	FD2K-99999-1	n/c
Frame Colors—Powder Coated		
<input type="radio"/> Black	<input type="radio"/> Pacific Blue	
<input type="radio"/> Blizzard Blue	<input type="radio"/> Passion Purple	
<input type="radio"/> Blue Blast	<input type="radio"/> Pink Parfait	
<input type="radio"/> Blue Stardust	<input type="radio"/> Really Red	
<input type="radio"/> Cosmic Purple	<input type="radio"/> Red Sunset	
<input type="radio"/> Extreme Green	<input type="radio"/> Sparkle Green	
<input type="radio"/> Mystic Teal	<input type="radio"/> Twilight	
<input type="radio"/> Octane Orange	<input type="radio"/> Yellow	
Frame Width		

Specify Width:

<input type="radio"/> <b>11" Wide</b>	FD2K-30210	n/c
<input type="radio"/> <b>12" Wide</b>	FD2K-30220	n/c
<input type="radio"/> <b>13" Wide</b>	FD2K-30230	n/c
<input type="radio"/> <b>14" Wide</b>	FD2K-30240	n/c
<input type="radio"/> <b>15" Wide</b>	FD2K-30290	n/c
<input type="radio"/> <b>16" Wide</b>	FD2K-30260	n/c
<input type="radio"/> <b>17" Wide</b>	FD2K-30330	n/c
<input type="radio"/> <b>18" Wide</b>	FD2K-30280	n/c

Width growth requires new cross braces and seat sling.

Frame Length		
<u>Specify Frame Length:</u>		
<input type="radio"/> <b>Regular (10"- 14")</b> 4" depth adjustment from 10" on regular frame. 3" additional growth achieved by replacing front frame section on regular frame only.	FD2K-30544	n/c
<b>Set-up Depth:</b>		
<input type="radio"/> <b>10" Deep</b>	FD2K-10	n/c
<input type="radio"/> <b>11" Deep</b>	FD2K-11	n/c
<input type="radio"/> <b>12" Deep</b>	FD2K-12	n/c
<input type="radio"/> <b>13" Deep</b>	FD2K-13	n/c
<input type="radio"/> <b>14" Deep</b>	FD2K-14	n/c
<input type="radio"/> <b>Long (13"- 17")</b> 4" depth adjustment from 13" on long frame.	FD2K-30545	<b>\$ 75</b>
<b>Set-up Depth:</b>		
<input type="radio"/> <b>13" Deep</b>	FD2K-13	n/c
<input type="radio"/> <b>14" Deep</b>	FD2K-14	n/c
<input type="radio"/> <b>15" Deep</b>	FD2K-15	n/c
<input type="radio"/> <b>16" Deep</b>	FD2K-16	n/c
<input type="radio"/> <b>17" Deep</b>	FD2K-17	n/c

Sling Seat		
<input type="radio"/> <b>Omit Seat Sling</b>	FD2K-39004	<b>(\$15)</b>
<input type="radio"/> <b>Sling Seat</b>	FD2K-39011	n/c
<u>Specify Sling Length:</u> Front of seat rail to back upholstery.		
<input type="radio"/> <b>10" Deep</b>	FD2K-30500	<b>\$ 55</b>
<input type="radio"/> <b>11" Deep</b>	FD2K-30501	<b>\$ 55</b>
<input type="radio"/> <b>12" Deep</b>	FD2K-30502	n/c
<input type="radio"/> <b>13" Deep</b>	FD2K-30503	<b>\$ 55</b>
<input type="radio"/> <b>14" Deep</b>	FD2K-30504	n/c
<input type="radio"/> <b>15" Deep</b>	FD2K-30505	<b>\$ 55</b>
<input type="radio"/> <b>16" Deep</b>	FD2K-30506	n/c
<input type="radio"/> <b>17" Deep</b>	FD2K-30507	<b>\$ 55</b>

Cushions		
<input type="radio"/> <b>1" Cushion</b>	FD2K-39003	<b>\$ 85</b>
<input type="radio"/> <b>2" Cushion</b>	FD2K-39007	<b>\$ 95</b>
Includes 4# foam & black knit cover.		

Date: \_\_\_\_\_ PO#: \_\_\_\_\_

Marked For \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

**Back Canes / Push Handles  
Hand Grip Bend Options**

- Straight** - No push handles  
With fold over sling back FD2K-30600 n/c
- Specify height from rear seat rail to top of straight cane:
- 10" High** FD2K-30610 n/c
- 11" High** FD2K-30611 n/c
- 12" High** FD2K-30612 n/c
- 13" High** FD2K-30613 n/c
- 14" High** FD2K-30614 n/c
- 15" High** FD2K-30615 n/c

- Traditional Bend Push Handles** FD2K-30518 n/c
- Ergo Bend Push Handles** FD2K-30519 n/c

- 8 Degree Back Cane Bend** FD2K-30511 \$ 25

Specify height from rear seat rail to top of bend on cane:

- Short** (19" High) FD2K-30515 n/c
- Medium** (23" High) FD2K-30514 n/c
- Tall** (26" High) FD2K-30516 n/c

**Back Angle**

Specify Degrees:

- 0 Degree** FD2K-39019 n/c
- 10 Degree** FD2K-39020 n/c

- Stroller Push Handle Extension - Quick Release** FD2K-30517 \$100

- One Piece Stroller Push Handle - Extension - Quick Release** FD2K-30526 \$ 125

- One Piece Push Button Angle Adjustable Stroller Handle** FD2K-30529 \$175

Back canes can achieve an additional height of 6.5" when quick release push handle extensions are ordered. All back cane heights are measured from seat rail to top of post.

**Sling Back**

- Omit Sling Back** FD2K-39005 n/c
- Standard Sling Back\*** FD2K-39006 n/c
- Tension Adjustable Back\*** FD2K-39008 \$231

Specify Sling Back Height:

Measured from top of seat to top of back upholstery.

- 10" High** FD2K-30910 n/c
- 11" High** FD2K-30911 n/c
- 12" High** FD2K-30912 n/c
- 13" High** FD2K-30913 n/c
- 14" High** FD2K-30914 n/c
- 15" High** FD2K-30915 n/c
- 16" High** FD2K-30916 n/c
- 17" High** FD2K-30917 n/c
- 18" High** FD2K-30918 n/c

\*Back upholstery height must be chosen when ordering sling or tension adjustable back upholstery

**Swing Away Hangers**

- Omit Hangers & Footplate** FD2K-33000 **(\$30)**
- Spring Loaded Pop Up** FD2K-30272 n/c
- Latch Style Swing Away** FD2K-30278 n/c

Specify Degrees:

- 90 Degree Hangers\*** FD2K-33334 n/c
- 75 Degree Hangers\*** FD2K-33335 n/c

\*See Configuration Chart



- Elevating Legrest** FD2K-30279 \$235

**Footrest Options**



- Omit Footplates Only** FD2K-36404 **(\$15)**
- Composite Footplate - Flip Up** FD2K-17070 n/c
- Kids Footplate 14w & Under Flip Up Footplate—Adjustable Angle** FD2K-36020 \$150  
For chairs 14" wide and under.
- Adult Footplate 15w & Over (E0951) Flip Up Footplate—Adjustable Angle w/ Heel Loop** FD2K-36021 \$150  
For chairs 15" wide and over.
- Omit Heel Loop** FD2K-36023 n/c
- Heel Loop (E0951)** FD2K-36022 \$ 50
- Footboard—Removable, Angle Adj.** FD2K-36029 \$125
- Hanger Extensions - .5" - 4.5" Long** FD2K-36402 \$ 55
- Hanger Extensions - 1.5" - 6.5" Long** FD2K-36403 \$ 65  
1.5" - 6.5" long only available with 75 degree hangers.

Configuration	Knee to Heel Range (from seat rail)
Hanger with High Mount Clamps	4" - 9"
Hanger Only	6" - 9"
Hanger with .5" - 4.5" Extension	to 13.5"
Hanger with 1.5" - 6.5" Extension	to 15.5"

- High Mount Clamp** FD2K-36407 \$100  
Not available on 75 degree hangers.

- Knee to Heel Measurement** \_\_\_\_\_ Inches

Date: \_\_\_\_\_ PO#: \_\_\_\_\_

Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

**Caster Wheels and Forks**

<input type="radio"/> Omit Caster Wheels	FD2K-34333	n/c
<input type="radio"/> 3" Poly Casters	FD2K-34888	\$ 55
<input type="radio"/> 3" Light Up Casters	FD2K-34111	\$ 55
3" casters recommended for low frame configuration.		
<input type="radio"/> 4" Poly Casters	FD2K-40000	\$ 65
<input type="radio"/> 4" Light Up Casters	FD2K-41000	\$ 99
<input type="radio"/> 5" Poly Casters	FD2K-34555	n/c
<input type="radio"/> 5" x 1 1/2" Poly Casters	FD2K-34222	\$ 85
<input type="radio"/> 6" Poly Casters	FD2K-34666	n/c
<input type="radio"/> 6" x 1.5" Poly Casters	FD2K-34999	\$ 85
<input type="radio"/> 6" x 2" Poly Casters with Tread	FD2K-34444	\$ 90
<input type="radio"/> 6" Pneumatic Casters	FD2K-34667	\$ 75
<input type="radio"/> 7" Poly Casters	FD2K-34777	\$ 75
<input type="radio"/> 7" x 1.5" Poly Casters	FD2K-34000	\$116
<input type="radio"/> 7" x 2" Poly Casters	FD2K-34200	\$116
<input type="radio"/> Frog Legs (E1015)	FD2K-35000	\$325
Not available with 6" x 2", 7", 7" x 1.5" & 7" x 2" casters.		
<input type="radio"/> Caster Stem Bolt - .75" Longer	FD2K-34778	n/c

**Anti-Tip Tubes**

<input type="radio"/> Anti-Tip Tubes—Fixed (E0971)	FD2K-31600	\$ 95
<input type="radio"/> Anti-Tip Tubes—Swing Away (E0971)	FD2K-31602	\$ 95
Can not be used with reverse configuration.		

**Wheels**

Patented double lock axle receiver included at no charge.

<input type="radio"/> Omit Wheels & Axles	FD2K-31524	(\$55)
---	------------	--------



Mag Wheel

**Mag Wheels**

<input type="radio"/> 12" Pneumatic Mag	FD2K-31525	\$ 40
Handrims are not available with 12" wheels.		

<input type="radio"/> 20" Pneumatic Mag	FD2K-31527	n/c
<input type="radio"/> 22" Pneumatic Mag	FD2K-31528	n/c
<input type="radio"/> 24" Pneumatic Mag	FD2K-31529	n/c

Aluminum handrims are included on 20", 22" and 24" wheels.  
To omit handrims, please check next column.  
24" wheels are not available with small frame lengths.

**Spoke Wheels**

<input type="radio"/> 20" Pneumatic Spoke	FD2K-31540	\$200
<input type="radio"/> 22" Pneumatic Spoke	FD2K-31541	\$200
<input type="radio"/> 24" Pneumatic Spoke	FD2K-31542	\$200

Aluminum handrims are included on 20", 22" and 24" wheels.  
To omit handrims, please check next column.  
24" wheels are not available with small frame lengths.  
Plastic coated handrims are not available with spoke wheels.  
Spoke wheels are not available with hub locks.



Spoke Wheel

**Wheel Accessories**

<input type="radio"/> Omit Handrims	FD2K-31006	(\$20)
<input type="radio"/> Aluminum Handrims 16" & 18" (pr.)	FD2K-31007	\$166
<b>Natural Fit Handrims</b> Please contact your customer service representative for details.		
<input type="radio"/> Plastic Coated Handrims (pr.)	FD2K-31005	\$117
Plastic coated handrims are not available with spoke wheels.		
<input type="radio"/> Airless Inserts (E2213)	FD2K-31000	\$ 75
Only available with pneumatic tire		
<input type="radio"/> Camber Plate	FD2K-31003	\$ 55
<u>Specify Degrees:</u>		
<input type="radio"/> 3 Degrees	FD2K-31011	n/c
<input type="radio"/> 6 Degrees	FD2K-31012	n/c

**Wheel Locks**

Choose type of lock - either standard or hub locks.

Specify lever type:

<input type="radio"/> Push to Lock*	FD2K-31501	n/c
<input type="radio"/> Pull to Lock*	FD2K-31506	n/c
<input type="radio"/> Rear Mount Push to Lock	FD2K-31507	n/c
<input type="radio"/> Rear Mount Pull to Lock	FD2K-31508	n/c
<input type="radio"/> Wheel Lock Handle Ext. (E0961)	FD2K-31801	\$ 60



\*If both wheel lock and hub lock are ordered, an "off chair" charge will be applied

OR

**Hub Locks**

(Not available with spoke wheels)

Specify lever location on frame:

Determines cable length and location of engagement lever for client access.

<input type="radio"/> Front Frame	FD2K-31509	n/c
<input type="radio"/> Middle Frame	FD2K-31510	n/c
<input type="radio"/> Rear Frame	FD2K-31511	n/c
<input type="radio"/> Back Canes	FD2K-31512	n/c

Date: \_\_\_\_\_ PO#: \_\_\_\_\_

Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

**Armrests**

- Omit Flip Back Armrest** FD2K-38209 **(\$20)**
- Back Cane Mount**
- Flip Back Armrest - 13" Long** FD2K-38204 **n/c**
- Flip Back Armrest - 16" Long** FD2K-38202 **n/c**
- Locking Flip Back Armrest - 13" Long (E0973)** FD2K-38201 **\$235**
- Locking Flip Back Armrest - 16" Long (E0973)** FD2K-38200 **\$235**
- Flip Back Armrests—Telescoping** FD2-38208 **\$100**
- Locking Flip Back Armrests—Telescoping** FD2-38207 **\$260**
- Frame Mount**
- Height Adjustable Padded Armrest (E0973)** FD2K-38203 **\$235**



Flip Back Armrest



Locking Flip Back Armrest



Height Adjustable Padded Armrest

- Desk Length Flat Armrest Pad** FD2K-38213 **n/c**
- Full Length Flat Armrest Pad** FD2K-38214 **n/c**
- Seat Rail to Top of Armrest** \_\_\_\_\_ **Inches**

**Accessories**

- Medical Necessities Bag (black only)** FD2K-23207 **\$ 67**
- Back Pack (black only)** FD2K-TFBP **\$ 55**
- Clothing Guards (fabric)** FD2K-CGF **\$120**
- Clothing Guards (Quick Release - rigid)** FD2K-CGR **\$120**
- Spoke Guards** FD2K-SG **\$ 90**

**\*Positioning Belts (E0978)**

- Streamline - Two Point Padded**
- Center Pull—Metal Push Button Buckle** SB-3332 **\$ 88**
- Center Pull—Metal Push Button with 1" Buckle & .75" Webbing** SB-3332-1 **\$ 88**
- Center Pull - EZ Release Push Button Buckle** SB-3334 **\$ 88**
- Streamline—Non-Padded**
- Center Pull—Metal Push Button Buckle** SB-3410 **\$ 29**
- Hipster - Two Point Padded**
- Center Pull—Metal Push Button Buckle** HIP-3002 **\$109**
- Center Pull - EZ Release Push Button Buckle** HIP-3004 **\$109**
- Hipster - Non-Padded**
- Center Pull—Metal Push Button** HIP-3021 **\$ 50**
- Center Pull—Metal Push Button w/ 1" Buckle & .75" Webbing** HIP-3021-1 **\$ 50**

\* Many other belts are available—the ones listed represent a few of our most popular. Reference our catalog or ask your customer service representative for more options.

**Wheel & Caster Matrix**

Front of Seat Rail To Floor	Rear of Seat Rail to Floor	Caster Size	Wheel Size	Fork Size	Arc Over Rail
15.5"	15.25"	3"	20"	Short	5.5"
15.5"	15.25"	3"	22"	Short	7.25"
17.25"	17"	3"	20"	Long	3.25"
17.25"	17"	3"	22"	Long	5.25"
17.5"	17.25"	5"	20"	Long	3.25"
17.5"	17.25"	5"	22"	Long	5.5"
17.5"	17.5"	5"	24"	Long	7"
18"	17.75"	6"	20"	Long	2.5"
18"	17.75"	6"	22"	Long	4.75"
18"	17.75"	6"	24"	Long	6.5"

**Maximum Rider Weight—250 lbs.**

Frame weight approximately 16.5 lbs.