



Freedom Designs, Inc.

FREEDOM 2

2241 N. Madera Rd. Simi Valley, CA 93065 Phone: (800) 331-8551 Fax (888) 582-1509 www.freedomdesigns.com

Bill to:	Date: _____ PO#: _____	Ship to:	Company: _____
	Company: _____		Address: _____
	Buyer: _____		City: _____ State: _____ Zip _____
	Marked For: _____		Phone: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Frame Type		
<input type="radio"/> Freedom 2 Depth Adjustable Rear Frame (E1236) or (K0005)	FD2-30270	\$2085
<u>Specify Frame Type:</u>		
<input type="radio"/> Folding	FD2-30275	n/c
<input type="radio"/> Rigid	FD2-30276	n/c
Transport Option		
<input type="radio"/> Transport Package Available with FDI WC19 approved seating or on frame only.	FD2-99998	\$150
<input type="radio"/> No Transport Package	FD2-99999-1	n/c
Frame Colors—Powder Coated		
<input type="radio"/> Black	<input type="radio"/> Pacific Blue	
<input type="radio"/> Blizzard Blue	<input type="radio"/> Passion Purple	
<input type="radio"/> Blue Blast	<input type="radio"/> Pink Parfait	
<input type="radio"/> Blue Stardust	<input type="radio"/> Really Red	
<input type="radio"/> Cosmic Purple	<input type="radio"/> Red Sunset	
<input type="radio"/> Extreme Green	<input type="radio"/> Sparkle Green	
<input type="radio"/> Mystic Teal	<input type="radio"/> Twilight	
<input type="radio"/> Octane Orange	<input type="radio"/> Yellow	
Frame Widths		
<u>Specify Width:</u>		
<input type="radio"/> 11" Wide	FD2-30210	n/c
<input type="radio"/> 12" Wide	FD2-30220	n/c
<input type="radio"/> 13" Wide	FD2-30230	n/c
<input type="radio"/> 14" Wide	FD2-30240	n/c
<input type="radio"/> 15" Wide	FD2-30290	n/c
<input type="radio"/> 16" Wide	FD2-30260	n/c
<input type="radio"/> 17" Wide	FD2-30330	n/c
<input type="radio"/> 18" Wide	FD2-30280	n/c
Width growth requires new cross braces and sling seat .		

Frame Length		
<u>Specify Frame Length:</u>		
<input type="radio"/> Regular (15"- 19") 4" depth adjustment from 15" on regular frame. 2" additional growth achieved by replacing front frame section on regular frame only.	FD2-30558	n/c
Set-up Depth:		
<input type="radio"/> 15" Deep	FD2-15	n/c
<input type="radio"/> 16" Deep	FD2-16	n/c
<input type="radio"/> 17" Deep	FD2-17	n/c
<input type="radio"/> 18" Deep	FD2-18	n/c
<input type="radio"/> 19" Deep	FD2-19	n/c
<input type="radio"/> Long (17"- 21") 4" depth adjustment from 17" on long frame.	FD2-30559	\$ 75
Set-up Depth:		
<input type="radio"/> 17" Deep	FD2-17	n/c
<input type="radio"/> 18" Deep	FD2-18	n/c
<input type="radio"/> 19" Deep	FD2-19	n/c
<input type="radio"/> 20" Deep	FD2-20	n/c
<input type="radio"/> 21" Deep	FD2-21	n/c
Sling Seat Length		
<input type="radio"/> Omit Seat Sling	FD2-39004	(\$15)
<input type="radio"/> Sling Seat	FD2-39011	n/c
<u>Specify Sling Length:</u> Front of seat rail to back upholstery.		
<input type="radio"/> 15" Deep	FD2-39015	\$ 55
<input type="radio"/> 16" Deep	FD2-39016	n/c
<input type="radio"/> 17" Deep	FD2-39017	\$ 55
<input type="radio"/> 18" Deep	FD2-39018	n/c
<input type="radio"/> 19" Deep	FD2-39021	\$ 75
<input type="radio"/> 20" Deep	FD2-39022	\$ 75
<input type="radio"/> 21" Deep	FD2-39023	\$ 75
Cushions		
<input type="radio"/> 1" Cushion	FD2-39003	\$ 85
<input type="radio"/> 2" Cushion Includes 4# foam & black knit cover.	FD2-39007	\$ 95



Date: _____ PO#: _____

Marked For _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

**Back Canes / Push Handles
Hand Grip Bend Options**

- Straight** - No push handles
With fold over sling back FD2-30600 n/c
- Specify height from rear seat rail to top of straight cane:
- 10" High** FD2-30610 n/c
 - 11" High** FD2-30611 n/c
 - 12" High** FD2-30612 n/c
 - 13" High** FD2-30613 n/c
 - 14" High** FD2-30614 n/c
 - 15" High** FD2-30615 n/c

- Traditional Bend Push Handles** FD2-30518 n/c
- Ergo Bend Push Handles** FD2-30519 n/c
- 8 Degree Back Cane Bend** FD2-30511 **\$ 25**

- Specify height from rear seat rail to top of bend on cane:
- Short (19" High)** FD2-30515 n/c
 - Medium (23" High)** FD2-30514 n/c
 - Tall (26" High)** FD2-30516 n/c

Back Angle

- Specify Degrees:
- 0 Degree** FD2-39019 n/c
 - 10 Degree** FD2-39020 n/c
 - Stroller Push Handle Extension - Quick Release** FD2-30517 **\$100**
 - One Piece Stroller Push Handle - Extension - Quick Release** FD2-30526 **\$125**
 - One Piece Push Button Angle Adjustable Stroller Handle** FD2-30529 **\$175**

Back canes can achieve an additional height of 6.5" when quick release push handle extensions are ordered. All back cane heights are measured from the seat rail to top of post.

Sling Back

- Omit Sling Back** FD2-39005 n/c
- Standard Sling Back*** FD2-39006 n/c
- Tension Adjustable Back*** FD2-39008 **\$231**

Specify Sling Back Height:
Measured from top of seat to top of back upholstery.

- 10" High** FD2-30910 n/c
- 11" High** FD2-30911 n/c
- 12" High** FD2-30912 n/c
- 13" High** FD2-30913 n/c
- 14" High** FD2-30914 n/c
- 15" High** FD2-30915 n/c
- 16" High** FD2-30916 n/c
- 17" High** FD2-30917 n/c
- 18" High** FD2-30918 n/c

*Back upholstery height must be chosen when ordering sling or tension adjustable back upholstery.

Footrest Hanger Options


- Omit Hangers & Footplate** FD2-33000 **(\$30)**
- Spring Loaded Pop Up** FD2-30272 n/c
- Latch Style Swing Away** FD2-30278 n/c

- Specify Degrees:
- 90 Degree Hangers*** FD2-33334 n/c
 - 75 Degree Hangers*** FD2-33335 n/c
- *See Configuration Chart



- Elevating Legrest (E0990)** FD2-30279 **\$235**
Not available with latch style.

Footrest Options



- Omit Footplates Only** FD2-36404 **(\$15)**
- Composite Footplate - Flip Up** FD2-17070 n/c
- Kids Footplate 14w & Under Flip Up Footplate—Adjustable Angle** FD2-36020 **\$150**
For chairs 14" wide and under.
- Adult Footplate 15w & Over (E0951) Flip Up Footplate—Adjustable Angle w/ Heel Loop** FD2-36021 **\$150**
For chairs 15" wide and over.
- Omit Heel Loop** FD2-36023 n/c
- Heel Loop (E0951)** FD2-36022 **\$ 50**
- Footboard—Removable, Angle Adj.** FD2-36029 **\$125**
- Hanger Extensions - .5" - 4.5" Long** FD2-36402 **\$ 55**
- Hanger Extensions - 1.5" - 6.5" Long** FD2-36403 **\$ 65**
1.5" - 6.5" long only available with 75 degree hangers.

Configuration	Knee to Heel Range (from seat rail)
Hanger with High Mount Clamps	4" - 9"
Hanger Only	6" - 9"
Hanger with .5" - 4.5" Extension	to 13 1/2"
Hanger with 1.5" - 6.5" Extension	to 15 1/2"


- High Mount Clamp** FD2-36407 **\$100**
Not available on 75 degree hangers.
- Knee to Heel Measurement** _____ **Inches**

Date: _____ PO#: _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Casters		
<input type="radio"/> Omit Caster Wheels	FD2-34333	n/c
<input type="radio"/> 3" Poly Casters	FD2-34888	\$ 55
<input type="radio"/> 3" Light Up Casters	FD2-34111	\$ 55
3" casters in combination with 20" wheels are <u>not</u> available. 3" casters recommended for low frame configuration.		
<input type="radio"/> 4" Poly Casters	FD2-40000	\$ 65
<input type="radio"/> 4" Light Up Casters	FD2-41000	\$ 99
<input type="radio"/> 5" Poly Casters	FD2-34555	n/c
<input type="radio"/> 5" x 1.5" Poly Casters	FD2-34222	\$ 85
<input type="radio"/> 6" Poly Casters	FD2-34666	n/c
<input type="radio"/> 6" x 1.5" Poly Casters	FD2-34999	\$ 85
<input type="radio"/> 6" x 2" Poly Casters with Tread	FD2-34444	\$ 90
<input type="radio"/> 6" Pneumatic Casters	FD2-34667	\$ 75
<input type="radio"/> 7" Poly Casters	FD2-34777	\$ 75
<input type="radio"/> 7" x 1.5" Poly Casters	FD2-34000	\$116
<input type="radio"/> 7" x 2" Poly Casters	FD2-34200	\$116
<input type="radio"/> Frog Legs (E1015)	FD2-35000	\$325
Not available with 6" x 2", 7", 7" x 1.5" & 7" x 2" casters.		
<input type="radio"/> Caster Stem Bolt - 75" Longer	FD2-34778	n/c
<input type="radio"/> Caster Stem Bolt - 1.5" Longer	FD2-34779	n/c

Anti-Tip Tubes		
<input type="radio"/> Anti-Tip Tubes—Fixed (E0971)	FD2-31600	\$ 95
<input type="radio"/> Anti-Tip Tubes—Swing Away (E0971)	FD2-31602	\$ 95
<u>Can not</u> be used with reverse configuration.		

Wheels		
Patented double lock axle receiver included at no charge.		
<input type="radio"/> Omit Wheels & Axles	FD2-31524	(\$55)
 Mag Wheels		
<input type="radio"/> 12" Pneumatic Mag	FD2-31525	\$ 40
Handrims <u>are not</u> available with 12" wheels.		
<input type="radio"/> 20" Pneumatic Mag	FD2-31527	n/c
<input type="radio"/> 22" Pneumatic Mag	FD2-31528	n/c
<input type="radio"/> 24" Pneumatic Mag	FD2-31529	n/c
Aluminum handrims are included on 20", 22" and 24" wheels. To <u>omit handrims</u> , please check next column. 24" wheels <u>are not</u> available with small frame lengths.		

Spoke Wheels		
		
<input type="radio"/> 20" Pneumatic Spoke	FD2-31540	\$200
<input type="radio"/> 22" Pneumatic Spoke	FD2-31541	\$200
<input type="radio"/> 24" Pneumatic Spoke	FD2-31542	\$200
Aluminum handrims are included on 20", 22" and 24" wheels. To <u>omit handrims</u> , please check next column. 24" wheels <u>are not</u> available with small frame lengths. Plastic coated handrims <u>are not</u> available with spoke wheels. Spoke wheels <u>are not</u> available with hub locks		

Wheel Accessories		
<input type="radio"/> Omit Handrims	FD2-31006	(\$20)
<input type="radio"/> Aluminum Handrims 16" & 18" (pr.)	FD2-31007	\$166
Natural Fit Handrims Please contact your customer service representative for details.		
<input type="radio"/> Plastic Coated Handrims (pr.)	FD2-31005	\$117
Plastic coated handrims <u>are not</u> available with spoke wheels.		
<input type="radio"/> Airless Inserts (E2213)	FD2-31000	\$ 75
Only available with pneumatic tire		
<input type="radio"/> Camber Plate	FD2-31003	\$ 55
<u>Specify Degrees:</u>		
<input type="radio"/> 3 Degrees	FD2-31011	n/c
<input type="radio"/> 6 Degrees	FD2-31012	n/c

Wheel Locks

Choose type of lock - either standard or hub locks.

		
<u>Specify Lever Option:</u>		
<input type="radio"/> Push to Lock*	FD2-31501	n/c
<input type="radio"/> Pull to Lock*	FD2-31506	n/c
<input type="radio"/> Rear Mount Push to Lock	FD2-31507	n/c
<input type="radio"/> Rear Mount Pull to Lock	FD2-31508	n/c
<input type="radio"/> Wheel Lock Handle Ext. (E0961)	FD2-31801	\$ 60

*If both wheel lock and hub lock are ordered, an "off chair" charge will be applied

OR

Hub Locks

(Not available with spoke wheels)

Specify lever location on frame
Determines cable length and location of engagement lever for client access.

<input type="radio"/> Front Frame	FD2-31509	n/c
<input type="radio"/> Middle Frame	FD2-31510	n/c
<input type="radio"/> Rear Frame	FD2-31511	n/c
<input type="radio"/> Back Canes	FD2-31512	n/c

Date: _____ PO#: _____

Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Armrests		
<input type="radio"/> Omit Flip Back Armrest	FD2-38209	(\$20)
Back Cane Mount		
<input type="radio"/> Flip Back Armrest - 13" Long	FD2-38204	n/c
<input type="radio"/> Flip Back Armrest - 16" Long	FD2-38202	n/c
<input type="radio"/> Locking Flip Back Armrest - 13" Long (E0973)	FD2-38201	\$235
<input type="radio"/> Locking Flip Back Armrest - 16" Long (E0973)	FD2-38200	\$235
<input type="radio"/> Flip Back Armrests—Telescoping	FD2-38208	\$100
<input type="radio"/> Locking Flip Back Armrests—Telescoping	FD2-38207	\$260
Frame Mount		
<input type="radio"/> Height Adjustable Padded Armrest (E0973)	FD2-38203	\$235
<input type="radio"/> Desk Length Flat Armrest Pad	FD2-38213	n/c
<input type="radio"/> Full Length Flat Armrest Pad Only available on 16" armrests.	FD2-38214	n/c
<input type="radio"/> Seat Rail to Top of Armrest _____ Inches		



Flip Back Armrest



Locking Flip Back Armrest



Height Adjustable Padded Armrest

*Pelvic Positioning Belts (E0978)		
Streamline - Two Point Padded		
<input type="radio"/> Center Pull—Metal Push Button Buckle	SB-3332	\$ 88
<input type="radio"/> Center Pull—Metal Push Button with 1" Buckle & .75" Webbing	SB-3332-1	\$ 88
<input type="radio"/> Center Pull - EZ Release Push Button Buckle	SB-3334	\$ 88
Streamline—Non-Padded		
<input type="radio"/> Center Pull—Metal Push Button Buckle	SB-3410	\$ 29
Hipster - Two Point Padded		
<input type="radio"/> Center Pull—Metal Push Button Buckle	HIP-3002	\$109
<input type="radio"/> Center Pull - EZ Release Push Button Buckle	HIP-3004	\$109
Hipster - Non-Padded		
<input type="radio"/> Center Pull—Metal Push Button	HIP-3021	\$ 50
<input type="radio"/> Center Pull—Metal Push Button w/ 1" Buckle & .75" Webbing	HIP-3021-1	\$ 50

* Many other belts are available—the ones listed represent a few of our most popular. Reference our catalog or ask your customer service representative for more options.

Wheel and Caster Matrix

Front of Seat Rail To Floor	Rear of Seat Rail to Floor	Caster Size	Wheel Size	Fork Size	Arc Over Rail
15.5"	15.25"	3"	20"	Short	5.5"
15.5"	15.25"	3"	22"	Short	7.25"
17.25"	17"	3"	20"	Long	3.25"
17.25"	17"	3"	22"	Long	5.25"
17.5"	17.25"	5"	20"	Long	3.25"
17.5"	17.25"	5"	22"	Long	5.5"
17.5"	17.5"	5"	24"	Long	7"
18"	17.75"	6"	20"	Long	2.5"
18"	17.75"	6"	22"	Long	4.75"
18"	17.75"	6"	24"	Long	6.5"

**Maximum Rider Weight—250 lbs.
Frame Weight Approximately 16.5 lbs.**

Accessories		
<input type="radio"/> Medical Necessities Bag (black only)	FD2-23207	\$ 67
<input type="radio"/> Back Pack (black only)	FD2-TFBP	\$ 55
<input type="radio"/> Clothing Guards (fabric)	FD2-CGF	\$120
<input type="radio"/> Clothing Guards (Quick Release - rigid)	FD2-CGR	\$120
<input type="radio"/> Spoke Guards	FD2-SG	\$ 90