



Freedom[®] NXT[™]

Generation NEXT Folding Tilt-in-Space Mobility

*Freedom: for Kids
(and other people too!)*



This year marks the 30th
Anniversary for
Freedom Designs, Inc.

To mark this occasion Freedom is
proud to introduce the NXT.
The NXT is the lightest, strongest
and easiest to adjust folding
tilt-in-space wheelchair yet.


The new Freedom NXT is truly the “Generation Next” in lightweight folding tilt-in-space mobility for Kids. You might even say it is the first “hybrid” by design.

The NXT incorporates all of the best qualities of rigid tilt frames like strength, durability and reliability with the premium qualities of folding tilt like lightweight, growth adjustment and compact portability.

Here are 30 reasons why the NXT is like no other folding tilt-in-space wheelchair in the world!



Our kids are our passion.™



- 1 Streamlined and aesthetically pleasing.
- 2 Compact to fit in small spaces.
- 3 First ever "rigid like" folding tilt.
- 4 Very light 22 lb. transfer weight.
- 5 Robust 200 lb. rider weight.
- 6 Easy to adjust and grow.
- 7 Eleven sizes in 1" increments.
- 8 Most sizes grow 2" in frame width.
- 9 Small, standard and large models.
- 10 Each model grows 4" in seat length.





11 Additional seat frame extensions in 2", 3" and 4" lengths.

12 Back angle adjust from 80° to 120°.

13 Shorter back post tube bracket.

14 Swing away or heavy duty pop-up foot rests.

15 Foot pedal activated tilt.

16 Single hand lever activated tilt.

17 Seat tube leveling guides.

18 Flared front frame.

19 Welded caster journals.

20 Aluminum mag front casters.

21

Lightweight front fork.

22

Universal rear axle plate.

23

Dual locking safety axle.

24

Transit tie downs mounted on base frame only.

25

Crash tested with Freedom Seating.

26

Combination anti-tipper and leverage bar.

27

Anti-rotation cross-brace tubes.

28

Anti-rotation seat frame connection.

29

Safety headlight reflecting coating.

30

Freedom Future Fit Program.



NXT™

Generation NEXT



Contact us:

Toll Free: (800) 331-8551

Local: (805) 582-0077

Fax: (888) 582-1509

Email:

CustomerService@freedomdesigns.com

2241 N. Madera Road

Simi Valley, CA 93065

www.freedomdesigns.com

FDI maintains a policy of continuous product improvement and reserves the right to change features, specifications and prices without prior notification.

Copyright © 2011 Freedom Designs, Inc.

All rights reserved. Printed in the USA.

